

HALL FLOWMETER KIT

10101.00.0.00

Key facts:

- The flowability of the powder can be ascertained quickly and precisely
- Different batches of powders can be compared
- A minimum of 10 measurements should be performed
- All components are safely stored in a handy suitcase

The Hall Flowmeter Kit (DIN EN ISO 4490) can ascertain the flowability of a powder quickly and precisely by determining the mobility of the powder particles in relation to each other. The measured value is expressed by the time required for 50 grams of metal powder to flow through a calibrated funnel.

The Hall Flowmeter Kit is stored in a suitcase that includes:

- Scale
- Stopwatch
- Measuring cup
- Comparative sand
- User manual

