

SCOTT VOLUMETER KIT

10102.00.0.00

Key facts:

- The test sample should have at least 100 cm³ in order to prepare a minimum of three individual samples
- The feed rate can be determined theoretically
- All components are safely stored in a handy suitcase

The Scott Volumeter Kit is designed to determine apparent density of the different powders as per DIN ISO 3923. Due to the different apparent density small deviations in particle size and structure can be detected.

The results of the measurements help to find out the powder feed rates.

The Scott Volumeter Kit is stored in a suitcase that also includes:

- Scale
- Measuring cup
- User manual

